

	Operating Mode	Réf :MO_Blot_Western Version 2
	<u>WESTERN BLOT</u>	date : 24 07 2018 page 1/2

1. Application field

Buffers for protein transfer and western-blotting

2. Referent persons

Véronique Santoni – Lionel Verdoucq

3. Health and safety :

Methanol :



[GHS02, GHS06, GHS08](#)

Hazard statements

[H225-H301 + H311 + H331-H370](#)

Personal Protective Equipment

[Eyeshields](#), [Faceshields](#), [full-face respirator \(US\)](#), [Gloves](#),

Thimerosal



[GHS06, GHS08, GHS09](#)

Hazard statements

[H300 + H310 + H330-H373-H410](#)

Personal Protective Equipment

[Eyeshields](#), [Faceshields](#), [full-face particle respirator type N100 \(US\)](#), [Gloves](#),

-Manipulate buffers with methanol under fume hood.

-Throw away effluents with methanol in specific waste « Effluents aqueux protéines + HPLC » .

4. Buffers :

Buffers for Protein transfer:

Depending on the aim and timing, the different buffers below can be used for protein blotting

1) Semi-dry blotting

glycine 39 mM 2,92 g

	Rédacteur	Vérificateur	Approbateur
Nom :	V Santoni L Verdoucq		C Tournaire

Tris 48 mM	5,81 g
SDS 0,0375 %	1,875 ml of 20 % solution
methanol 20 %	200 ml
H ₂ O qsp	1 l

*transfer at 0,8 mM/cm² for 1h 15 or
transfer at 12 V for 15 min , and 24 V for 45 min*

2) liquid transfer in CAPS-methanol

Caps 10 mM
methanol 10 % 100 ml
adjust to pH 11 with NaOH 6N
transfer 35 V for 16 h at 4°C

3) Liquid Tris-Glycin-methanol

Tris 10 mM 1.82 g /1.5l
Glycin 100 mM 11.26 g/1.5l
Methanol 10 % 150 ml
H₂O qsp 1.5 l
Transfer 70 v for 1h30 at 4°C

4) Liquid Tris-borate (for sequencing)

- equilibrate gels for 2 h in boric acid 50 mM, SDS 0,1 %
Adjust to pH 8 with NaOH 6N
- transfer in boric acid 50 mM
Tris 50 mM
35 V for 16 h at 4°C

Buffers for Western

PBS 10X buffer

KH ₂ PO ₄	40 mM
Na ₂ HPO ₄	160 mM
NaCl	1.15 M

Adjust pH to 7,3 with NaOH 6N

Add thimerosal 0,002 % (p/v) for long conservation

	Rédacteur	Vérificateur	Approbateur
Nom :	V Santoni L Verdoucq		C Tournaire

	Operating Mode	Réf :MO_Blot_Western Version 2
	<u>WESTERN BLOT</u>	date : 24 07 2018 page 3/2

PBSTB buffer

	finale conc	quantity
PBS	1X	50 ml of 10 X buffer
Tween-20	0,1 %	500 mg
BSA	1 %	5 g
H2O qsp		500 ml

Buffer volumes : 40 ml for small membrane in square petri dish and 80 ml for big membrane;

Possible to incubate primary Antibody in a sealed bag.

5. Operating mode for Western-blot :

After blotting, wash the PVDF membrane in PBS

Membrane saturation :

- In PBSTB

1 h at room temperature or ON at 4°C

- Primary antibody

Dilute in PBSTB, incubate 2 h at room temperature or ON at 4°C

- Wash with PBSTB 2 x 10 min

- Secondary antibody: peroxidase conjugate

dilute 20.000 X in PBSTB, incubate 1 h at room temperature

- Wash with PBS 4 x 5 min

- revelation: chemiluminescence

Drain the membrane on a paper

Incubate the membrane in presence of a mix of substrates "Super signal West Pico Chemiluminescent substrate" (Thermo Sci.):

x ml solution A + x ml solution B during 5 min

Drain the membrane on a paper

	Rédacteur	Vérificateur	Approbateur
Nom :	V Santoni L Verdoucq		C Tournaire

	Operating Mode	Réf :MO_Blot_Western Version 2
	<u>WESTERN BLOT</u>	date : 24 07 2018 page 4/2

Expose and reveal using LAS 3000.

	Rédacteur	Vérificateur	Approbateur
Nom :	V Santoni L Verdoucq		C Tournaire

	Operating Mode	Réf :MO_Blot_Western Version 2 date : 24 07 2018 page 5/2
	<u>WESTERN BLOT</u>	